

## Method Calibration

Name	Value
Section: Method	
Category: Identification	
ID	PH CAL
Type	Calibration
Version	5
Internal ID	M22
Name	Calibration
Section: Sample (Calib) (C1)	
Category: Sensor	
Sensor type	pH
Sensor	DG115-SC
Unit	pH
Category: Calibration	
Action	Calibration
Automatic buffer recognition	No
pH buffer list	METTLER TOLEDO USA (Ref. 25°C)
Calibration	Linear
Number of buffers	3
Buffer 1	4.01
Buffer 2	7.00
Buffer 3	10.01
Category: Temperature	
Temperature [°C]	25.0
Section: Titration stand (Titrationstand1)	
Category: Titration Stand	
Type	InMotion T/Tower A
Titration stand	InMotion T/1A
Head position	Sample
Lid handling	No
Section: Stir (Stir1)	
Category: Stir	
Speed [%]	30
Duration [s]	60
Condition	No
Section: Measure (normal) [1] (Measurenormal1)	
Category: Sensor	
Type	pH
Sensor	DG115-SC
Unit	mV
Category: Temperature Acquisition	
Temperature Acquisition	No
Category: Stir	

Name	Value
Speed [%]	30
Category: Acquisition of Measured Values	
Acquisition	Equilibrium controlled
dE [mV]	0.5
dt [s]	1
t(min) [s]	10
t(max) [s]	30
Mean value	No
Category: Condition	
Condition	No
Section: Record (Record1)	
Category: Report	
Report template	Calibration report
Print	No
Condition	No
Section: Calibration (Calibration1)	
Category: Calibration	
Sensor	DG115-SC
Limits	No
Condition	No
Section: Park (Park1)	
Category: Park	
Titration stand	InMotion T/1A
Position	Conditioning beaker
Condition	No
Section: Record (Record2)	
Category: Report	
Report template	Calibration report
Print	No
Condition	No





