

Method Silbiger TA (EP)r1

Name	Value
Section: Method	
Category: Identification	
ID	00002r1
Type	General Titration
Version	8
Internal ID	M55
Name	Silbiger TA (EP)r1
Section: Sample (S1)	
Category: Sample	
Sample type	Sample
Number of IDs	1
ID 1	--
Entry type	Weight
Lower limit [g]	59
Upper limit [g]	61
Density [g/mL]	1.0
Correction factor	0
Temperature [°C]	25.0
Entry	Arbitrary
InMotion reader	None
Number of sample factors	0
Section: Titration stand (Titrationstand1)	
Category: Titration Stand	
Type	InMotion T/Tower A
Titration stand	InMotion T/1A
Head position	Sample
Lid handling	No
Section: Stir (Stir1)	
Category: Stir	
Speed [%]	70
Duration [s]	60
Condition	No
Section: Titration (EP) [1] (TitrationEP1)	
Category: Titrant	
Titrant	HCl
Concentration [mol/L]	0.1
Category: Sensor	
Type	pH
Sensor	DG115-SC
Unit	mV
Category: Temperature Acquisition	
Temperature Acquisition	Yes

Name		Value
	Temperature sensor	DT1000
	Temperature unit	°C
Category: Stir		
	Speed [%]	30
Category: Predispose		
	Mode	Potential
	Potential [mV]	185
	Wait time [s]	240
Category: Control		
	Endpoint type	Absolute
	Tendency	None
	Endpoint value [mV]	235.0
	Control band [mV]	50
	Dosing rate (max) [mL/min]	0.5
	Dosing rate (min) [µL/min]	10
Category: Termination		
	At EP	Yes
	Termination delay [s]	0
	At Vmax [mL]	4.0
	Max. time infinite	Yes
	Max. time [s]	∞
Category: Accompanying Stating		
	Accompanying stating	No
Category: Condition		
	Condition	No
Section: Content (R3)		
Category: Calculation		
	Result	Content
	Result unit	mL
	Formula	R3=VEQ
	Constant C=	1
	M	M[Hydrochloric acid]
	z	z[Hydrochloric acid]
	Decimal places	3
	Result limits	No
	Extra statistical functions	No
	Send to buffer	No
	Write to Smart Tag	None
	Condition	No
Section: Content (R1)		
Category: Calculation		
	Result	Content
	Result unit	g

Name	Value
Formula	$R1=m$
Constant C=	$M/(10 \cdot z)$
M	M[Hydrochloric acid]
z	z[Hydrochloric acid]
Decimal places	3
Result limits	No
Extra statistical functions	No
Send to buffer	No
Write to Smart Tag	None
Condition	No

Section: Record (Record1)

Category: Report

Report template	Titration report
Print	No
Condition	No



